

Analytical Validation of the Thermo Fisher ABI HIV-1 Genotyping Kit for Illumina Sequencing, Benchmarked Against Sanger, Using the Exatype Platform

Executive Summary

[Hyrax Biosciences](#), through its [Exatype platform](#), delivers scalable HIV drug-resistance analysis to global programmes - processing over **150,000 HIV sequences** for **1,146 users** across **70+ organisations** in **103 countries**. More than half of these programmes operate in Global Health Equity settings, positioning [Exatype](#) as one of the most **widely adopted** and **trusted platforms** for routine Sanger-based HIV genotyping worldwide.

Building on this foundation, [Hyrax Biosciences](#) validated whether Illumina NGS data generated using the Thermo Fisher HIV-1 Genotyping Kit with Integrase, when analysed on [Exatype](#), delivers drug-resistance results equivalent to the Sanger (capillary electrophoresis sequencing) gold standard. This analysis shows **strong cross-platform concordance**, demonstrating that the kit performs consistently whether sequenced on Illumina or Sanger. [Exatype](#) was used as the **common analytical pipeline**, ensuring **standardised processing** and **fully comparable outputs**.

Key Outcomes:

- **Specificity:** >98.7% concordance between Illumina and Sanger consensus sequences.
 - **Repeatability:** ≥99.0% intra-sample concordance across replicates.
 - **Sensitivity:** Equivalent detection of mutations at ≥20% prevalence (Sanger gold standard); Illumina additionally resolves variants below this threshold.
 - **Throughput:** Illumina generated ~340,000 reads per replicate, delivering high coverage across PR, RT and IN targets.
-

Background and Rationale

At the end of 2022, **29.8 million people** were on antiretroviral therapy (World Health Organisation, HIV Statistics, 2025)¹, yet **HIV drug resistance** remains a major global health threat. Global HIV drug-resistance monitoring requires **validated workflows** that are **accurate, scalable, and compatible across sequencing technologies**. While Sanger sequencing

remains widely used, laboratories are increasingly shifting toward NGS platforms such as Illumina to improve throughput and sensitivity.

[Hyra Biosciences](#) undertook this **analytical validation** to assess whether Illumina sequencing data generated using the Thermo Fisher ABI HIV-1 Genotyping Kit with Integrase, analysed on the [Exatype platform](#), provides results consistent with Sanger sequencing while enabling deeper variant profiling. Both the Sanger sequencing and Illumina NGS data were processed using [Exatype](#) ensuring directly comparable outputs.

Validation Results

Study design: included 5 samples × 4 replicates for PR (protease) and RT (reverse transcriptase) - sequenced as a single amplicon - and IN (integrase). All samples and replicates were amplified with the Thermo Fisher ABI HIV-1 Genotyping Kit with Integrase, sequenced on Illumina and Sanger, and analysed on [Exatype](#).

Quantitative Outcomes




Figure 1: Intra- and Between-Sample Consensus Sequence Concordance

Metric	Result	Benchmark / Notes
Specificity	>98.7% Illumina-Sanger consensus concordance (Fig. A)	Confirms platform equivalency
Repeatability	>99.0% concordance across replicates (Fig. B)	Demonstrates high robustness
Sensitivity	Matching detection of major variants (≥20% Sanger standard)	Equivalent performance to Sanger
Sensitivity threshold	Sanger threshold 15 - 20% ² ; Illumina threshold ≥1%	Enhanced minor-variant resolution with Illumina

Read depth	~340,000 reads per replicate; ≥89% mapped	Supports high-confidence variant calls
Coverage	>15,000× average across PR, RT, IN	Exceeds minimum required (≥14x) for DR calling
Drug-resistance calls	Consistent across platforms	Reliable HIV-DR detection

Key observations:

- **Reproducibility** across Illumina and Sanger sequences was ≥98%, confirming **assay performance**.
- **Minor differences** between Illumina and Sanger (and among replicates) were primarily attributable to shifts in variant prevalence at **mixture sites**.

Analytical and Applied Implications enabled by Hyrax Biosciences' Exatype Platform

This **validation** confirms that programmes can continue using the same Thermo Fisher ABI HIV-1 Genotyping Kit they have relied on for years, while seamlessly expanding from Sanger to Illumina without changing their **analytical workflow**. By processing both Sanger and Illumina data through the same [Exatype](#) pipeline, laboratories retain identical interpretation rules, mixture thresholds, reporting formats, and drug-resistance outputs - **ensuring continuity across longitudinal datasets**.

Illumina's high read depth, combined with [Exatype's](#) automated QC and unified variant-calling framework, offers improved minor-variant detection without altering established Sanger-based thresholds. This **cross-platform consistency** allows laboratories to increase throughput and sensitivity using familiar reagents and a familiar analysis environment, **removing the need for new bioinformatics tools or retraining**. In practice, it reduces operational burden, supports scalable national surveillance, and keeps historical Sanger data comparable as programmes transition to NGS.

Conclusion

Illumina sequencing, generated using the Thermo Fisher ABI HIV-1 Genotyping Kit with Integrase and analysed through [Hyrax Biosciences' Exatype platform](#), provides **HIV-1 drug-resistance** results equivalent to Sanger while adding the benefits of high throughput and improved minor-variant resolution. This **analytical validation** confirms [Exatype](#) as a **robust, cross-platform solution** capable of supporting **large-scale HIV-DR surveillance and laboratory implementation across diverse sequencing environments**.



References

1. [World Health Organisation \(WHO\). HIV statistics, globally and by WHO region, Information Sheet, 2025](#)
2. [Steegan et al., 2023, Advancing HIV Drug Resistance Technologies and Strategies: Insights from South Africa's Experience and Future Directions for Resource-Limited Settings. Diagnostics, 13\(13\):2209](#)

